



Variable
Crimpable
Hook



Acc Kit



Square
Thread

Abso Chain



Dentos Inc.

For better orthodontic solutions

41, Seongseo-ro 68-gil, Dalseo-gu, Daegu, Korea 704-900
telephone. +82-(0)53-424-5947, +82-(0)53-592-5908
fax. +82-(0)53-592-5909

E-mail. webmaster@dentos.co.kr
www.dentos.co.kr

Dentos



Square
Thread



Crimpable
Stop

Variable
Crimpable
Hook



Coil Spring



























Abso Chain



Orthodontic Accessories

Accessories
- 5th ed (May. 2014)

NT(Nickel-Titanium) Closed Coil Spring

Length	Hole Diameter Ø	Strength			
		Ultra Light(007x036)	Light (009x036)	Medium (010x036)	Heavy(010x036)
8mm 	Ø 1.5	 NT15-8U	 NT15-8L	 NT15-8M	 NT15-8H
	Ø 2.0	 NT20-8U	 NT20-8L	 NT20-8M	 NT20-8H
	Ø 2.5	 NT25-8U	 NT25-8L	 NT25-8M	 NT25-8H
13mm 	Ø 1.5	 NT15-13U	 NT15-13L	 NT15-13M	 NT15-13H
	Ø 2.0	 NT20-13U	 NT20-13L	 NT20-13M	 NT20-13H
	Ø 2.5	 NT25-13U	 NT25-13L	 NT25-13M	 NT25-13H

Wire diameter

- Ultra Light 0.19mm (0.007")
- Light 0.21mm (0.009")
- Medium & Heavy 0.23mm (0.010")

Spring inner diameter

- 0.9mm (0.036")

Strength

- Ultra Light (U), Light (L), Medium (M), Heavy (H)


Length

- 8mm & 13mm

Hole diameter


- Ø1.5, Ø2.0, Ø2.5

Use with MI series




- NT 15 : use with SH, LH, PH series, GH 13~16 series, GH 20 series
- NT 20 : use with CH, DH series
- NT 25 : use with BH, TH, FH, HH series, GH 17~18 series
- 1 box contains 10 springs

Weight Range (8mm length)



Weight Range (13mm length)





NT Closed Coil Spring Case



◀ Each case is detachable




Coil Spring Holder (CSH)


















NT(Nickel-Titanium) Open Coil Spring

Wire sizes available : 022



Abso Chain

Type	Clear Type		Gray Type	
	0.4	0.5	0.4	0.5
 <p>Narrow</p>	 <p>N4-C</p> <p>0.4mm</p>	 <p>N5-C</p> <p>0.5mm</p>	 <p>N4-G</p> <p>0.4mm</p>	 <p>N5-G</p> <p>0.5mm</p>
 <p>Medium</p>	 <p>M4-C</p> <p>3mm</p> <p>0.4mm</p>	 <p>M5-C</p> <p>0.5mm</p>	 <p>M4-G</p> <p>0.4mm</p>	 <p>M5-G</p> <p>0.5mm</p>
 <p>Wide</p>	 <p>W4-C</p> <p>5mm</p> <p>0.4mm</p>	 <p>W5-C</p> <p>0.5mm</p>	 <p>W4-G</p> <p>0.4mm</p>	 <p>W5-G</p> <p>0.5mm</p>





Abso Chain



Square Thread

Square Thread







	A	B	C
ST-055	0.55mm	0.55mm	0.74mm
ST-045	0.45mm	0.45mm	0.67mm

Hook

Variable Crimpable Hooks

Wire sizes available : 016 ~ 022 x 028





"V" notch prevent rotation of hook during grip


- Grip the hook with the open part of the crimp facing externally by use of plier
- Insert the open part of crimp into orthodontic archwire while moving plier
- Apply force to the plier for squeezing the crimpable part
- Remove the plier from the crimpable hook
- You can cut into the proper length of hook with a cutter.

Crimpable Stop (VCH-CS)


Wire sizes available : 016 ~ 022 x 028



Crimpable Hook Plier (CHP)



Hook Case (H-Case)



Acc Kit

Components

1. NT Closed Coil Spring : 4pkg
(NT15-8M, NT15-13M, NT25-8M, NT25-13M)
2. Hook : 4pkg
(VCH-L-2, VCH-R-2, VCH-L-3, VCH-R-3)
3. Square Thread : 2ea
(ST-045, ST-055)
4. Abso Chain : 4ea
(N4-C, N5-C, M4-C, M5-C)
5. CHP : 1ea

Free

6. CSH : 1ea
7. H-Case : 1ea



Dentos Distributors

Japan

Shofu Inc.
33rd floor 3-16-2 Yushima, Bunkyo-ku, Tokyo [113-0034], JAPAN
Tel : +81-3-3832-1824 Fax : +81-3-3832-7682
info@shofu.co.jp

Taiwan

Shin Yean Trading Co., Ltd.
No 2, Alley 2, Lane 147, Sec.1, Nan Chang RD, Taipei, TAIWAN
Tel : +886-2-2391-8481 Fax : +886-2-2321-6986
shinyean@ms18.hinet.net

China

BEIJING FOCUS CORPORATION.
6F Fairyland Plaza, 83# Fuxing Rd, Haidian District, Beijing
100039, P.R. CHINA
Tel : +86-10-8861-2661 Fax : +86-10-8861-2690
wanglei@bjfic.com

Thailand

Dentapex Co., Ltd.
128/146 Payatai Plaza BD 13F, Payatai Road, Rajthevee,
Bangkok, 10400, THAILAND
Tel : +66-(0)2-216-1688 Fax : +66-(0)2-612-9042
walaiporn.g@dentapex.co.th

HongKong

Forever Green Dental Products Ltd.
Rm 1301, Manulife Provident Funds Place, 345 Nathan Road,
Kowloon, HONG KONG
Tel : +852-2388-2798 Fax : +852-2332-8183
forevergreendental@gmail.com

India / UAE

Dentos India Pvt. Ltd.
301, Lavlesh Court, Pandit Varde Road, Bandra (West) Mumbai,
400 050 INDIA
Tel : +91-22-4250-5100 Fax : +91-22-4250-5101
mail@dentosindia.com

Singapore / Indonesia

Fondaco Pte. Ltd.
7 Kaki Bukit Road 1, #03-06, Eunos Technolink 415937
SINGAPORE
Tel : +65-6392-2806 Fax : +65-6392-1296
fondaco@fondaco.sg

Malaysia

Denta Medic Technology
No. 38A, Jalan Renang 13/26, Section 13, 40100 Shah Alam,
Selangor, MALAYSIA.
Tel : +60-3-5510-4346 Fax : +60-3-5510-4345
info@denta-medic.com

Philippines

Ahescor Inc.
No. 5029 Filmore St., Palanan, 1235, Makati City, Manila,
PHILIPPINES
Tel : +63-2-551-4679 Fax : +63-2-551-4747
ahescor@gmail.com

Vietnam

D.O.E (Dentistry Of Excellence)
199 Nguyen Thi Minh Khai, Dist 1, Ho Chi Minh City, VIETNAM
Tel : +84-8-3925-3703 Fax : +84-8-3925-3700
doevietnam@gmail.com

Azerbaijan

RealDent LLC.
Senyushkin Str. 13, AZ 1025, Baku, AZERBAIJAN Republic
Tel : +99-412-490-6845 Fax : +99-412-490-2389
info@realdent.net

Pakistan / Afghanistan

SPIRAL INTERNATIONAL
Flat # 12-A, 2nd Floor, HBL Bldg, Agro flats, Shadman Market,
Lahore, 54000, PAKISTAN
Tel : +92-42-7597320 Fax : +92-42-7503699
spiralkp@yahoo.com

UK

Torque Orthodontics Ltd.
The Tannery, Home Farm, Esholt, Bradford, West Yorkshire, BD17
7RH, U.K
Tel : +44-1274-581058
john@torqueorthodontics.com

Spain

IMPER-ORTHO
Plaza Mayor, 40 - 2º - 2ª, 46600-Alzira (Valencia), SPAIN
Tel : +34-962-415-647 / 962-404-600
Fax : +34-962-401-289
info@imperortho.com

France / Israel / Belgium / Tunisia / Morocco

DENTAL FORCE
32, bd Paul Vaillant-Couturier, F-93100 MONTREUIL, FRANCE
Tel : +33-(0)1-4360-7708 Fax : +33-(0)1-4360-7905
info@dentalforce.fr

Netherlands

Dentsply Lomberg B.V.
Vredehofstraat 3-5, 3761 HA SOEST, [Post box 164],
NETHERLANDS
Tel : +31-(0)35-609-4753 Fax : +31-(0)35-602-5480
lomberg.order@dentsply.com

Swiss / Liechtenstein

FEANRO Ltd.
Gottfried Keller-Strasse 7, CH-8001, Zurich, SWITZERLAND
Tel : +41-44-562-1862 Fax : +41-44-562-1861
info@feanro.com

Poland

ORTO TRADING
ul. Husarii 11 B, 02-951, Warszawa, POLAND
Tel : +48-22-651-8409 Fax : +48-22-842-5368
ortotrading@ortotrading.com.pl

Turkey

Mediko Dental
Tevsikiye Cad. No:47 3/3, 34365, Nisantasi, Istanbul, TURKEY
Tel : +90-212-219-9755 Fax : +90-212-219-9760
info@medikodental.com

Greece / Cyprus

MULTIADIS VITSARPOPOULOS S.A.
348 MESOGION AVENUE, 153 41 AGIA, PARASKEVI, ATHENS,
GREECE
Tel : +30-210-65-41-340 Fax : +30-210-65-41-618
info@vitsaropoulos.gr

Hungary

Dr. Volom Dental & Medical Kft
1055 Budapest, Nyugati ter 9. I./3. Hungary
Tel : +36-1-311-6584 Fax : +36-1-489-3709
shop@drvolomdental.hu

Bulgaria

DENTALAB.BG Ltd.
33, Gen. Totleben Blvd. Sofia, 1606 Bulgaria
Tel : +359-2-423-4856 / +359-87-8725185
Fax : +359-2-423-4856 office@dentallab.bg

Argentina / Chile

PRODENTOS S.A
Terrero 1960 5º piso A, C.P. 1416 BMR, Capital Federal Buenos
Aires, ARGENTINA
Tel : +54-11-4139-4274 (4139-4295) mhprodent@yahoo.com.ar

Colombia

DENTOS COLOMBIA SAS
CRA 17 # 108 - 29 AP 100, BOGOTA - COLOMBIA
Tel : +57-(1)-4953673 / +57 (300) 4873039
juanderosabalen@gmail.com

Uruguay

Pro3implant S.R.L.
Colonia 2099, Floor 1st C.P. 11200, Montevideo, URUGUAY
Tel : +598-2408-3003 Fax : +598-2403-4163
info@pro3implant.com.uy

Panama

ImpOrt Corp.
Plaza Mont Blanc, Local #5, Avenida Abel Bravo Y Calle 60
Urbanizacion Obarrio, Panama, PANAMA
Tel : +507-263-0200 Fax : +507-264-5188
rm@frenos.cc

Dominican Rep.

Selladores Cobian S.A.
Gustavo Mejia Ricart #37, Edif. Boyero III, Apto. 202, Ensanche
Naco, Santo Domingo, República DOMINICANA
Tel : +1-809-547-1220 Fax : +1-809-549-3595
franciscmartino@gmail.com

Peru / Bolivia

SERVICIOS ODONTOLÓGICOS DEL PLATA E.I.R.L.
Ave. DE LAS ARTES 1482, San Borja, Lima 41, PERU
Tel : +51-1-436-4432 rafaelberenguel@gmail.com

Venezuela

L.V.M Dental Product C.A
Calle Rio Caroni, edificio El Moruco, piso 4,
apartamento 4-B. Cumbres de Curumo, Caracas 1080,
VENEZUELA
Tel : +58-414-030-0537 Fax : +58-414-125-0197
lucyperez30@yahoo.com

USA

Great Lakes Orthodontics
200 Cooper Ave. Tonawanda, N.Y 14150, USA
Tel : +1-800-828-7626 / 866-781-5075 Fax : +1-716-871-0550
joepanara@greatlakesortho.com

Australia

Dentos Australia
3/14 ETHEL ST. EASTWOOD, NSW 2122, AUSTRALIA
Tel : +61-(0)2-9874-6117 Mob.: 0416-028-014
Fax : +61-(0)2-9874-8922
dentos1@optusnet.com.au

Russia

NORTA Co. Ltd.
212, entr. № 5, 0A, Dibunovskaya str., Saint-Petersburg,
197183, RUSSIA
Tel : +7-812-935-8448 Fax : +7-812-640-2840
info@nortonline.com

Ukraine

Ortho-Medina
1013, 172 Gorkogo str. 03680 KIEV, UKRAINE
Tel : +38-(0)44-219-1503 Fax : +38-(0)44-219-1504
rezai_hamid@yahoo.com

Kazakhstan / Uzbekistan

NUR RAFFA Co., Ltd.
Altyn Gasyr Business Center office 112, Shevchenko street, 118,
Almaty, KAZAKHSTAN
Tel : +7-727-244-0459
uzshin@hanmail.net

Lebanon

Star Science Lebanon S.A.L
Gema Center, 3rd Floor, Dbayeh highway, Beirut, LEBANON
Tel : +961-3-10-1211 Fax : +961-4-54-3496
starsciencelb@gmail.com

Syria

Al Lahham Corp.
P.O. box 3727, karameh street, Abou Rummaneh, Damascus,
SYRIA
Tel : +963-11-33-33719 Fax : +963-11-33-13719
laham.cntr@yahoo.com

Kuwait

Dental Tech
Block# 9, Street# 3, Unit#16, Al-Jabriya, KUWAIT
Tel : +965-901-0185 Fax : +965-787-8704
leosin72@hotmail.com

Jordan

Al-Eiman Drug Store
7th circle, Airport st. beside city terminal, Amman, 11810.
JORDAN
Tel : +962-6-581-0108 Fax : +962-6-581-8850
samerissa@eiman-eds.com

Bahrain

Transalign Trading
Building 1135, Flat NO:11, STROAD: 3433, Block:734,
AALII TOWN,
[P.O Box 1282 - Manama] KINGDOM OF BAHRAIN
Tel : +973 13 631631 Fax : +973 17 911599
info@transalign.net

Oman

AL FARSI MEDICAL SUPPLIES
Office No.4, Ground Floor, Al Khuwair Str. Muscat, Boshier, Postal
Code 133, (P.O. BOX 158) Sultanate of OMAN
Tel : +968-2448-5625 Fax : +968-2448-5675
sales@alfarsi.me

Iraq

ECHO For Dental Equipments
Hospital Roundabout, Mosul, IRAQ
Tel : +964 7481752665 (Mob. 7707488499)
echodent@echodent.com

Egypt

ELMO Medical
511 Haram St., Nasr building, Giza Square, Giza, EGYPT
Tel : +20-1001031119 / +20-1111801107
Fax : +20-643345325
elmotayammed@gmail.com